

HEAD LUX

FIBEROPTIC & LED OPTIONS

The extended family of Ecleris Headlights,
Everlux & Headlux Products

Premium quality dual adjustment head band
Optimized illumination



FIBEROPTIC HEADLUX

LED HEADLUX



HEAD LUX

FIBEROPTIC & LED OPTIONS

The extended family of Ecleris Headlights,
Everlux & Headlux Products



High light transmission
Adjustable focus with extremely precise beam
Light weight – adjustable – comfortable

• FIBEROPTIC HEADLUX



Standard fiber optic
cable connection
(Storz type)

- Cold light
- Best combination with Ecleris Ever Lux Mini



• LED HEADLUX



Battery operated:
10 hours continuous
use at maximum
intensity

- LED life > 50.000 hours
- Rechargeable battery
- Color temperature 6.000 ° K
- Charger included



www.ecleris.com

ECLERIS
INNOVATION VALUE